

Dogic Lab Protocol for Purifying Kinesin Proteins

This protocol is for purifying His-tagged, truncated kinesin constructs using an AKTA.

Buffers

Wash Buffer:

50mM PIPES
4mM MgCl₂
50uM ATP
10mM β ME (4°)
20mM Imidazole

pH with KOH to 7.2

Elution Buffer:

50mM PIPES
4mM MgCl₂
50uM ATP
2mM DTT (4°)
500mM Imidazole

pH with HCl to 7.2

Lysing

- Remove cells from freezer and put on ice
- add protease inhibitors to wash buffer to make lysis buffer
 - 1 tablet/10mL of Pierce Protease Inhibitor Mini Tablet (Thermo #88665)
 - a 1mg/mL of PMSF (MP Biomedicals, # 195381)
 - add 1mg/mL of Lysozyme, Chicken, RZ3 Muramidase (US Biological # L9200)
- add 1-1.5x pellet volume of lysis buffer
- leave on ice to thaw
- resuspend and
- lyse cells using fine tip: 50% power, 50 % time, 15 minutes (Omni International Sonic Ruptor)
- every 5 minutes go in to re-support the falcon tube with ice (the lyser generates heat which melts the ice and can damage the proteins).
- centrifuge at $\geq 100,000$ Gs for ≥ 30 min at 4° to clarify lysate

Akta

- open software
- set pressure alarm (listed under column information)
- start flowing ethanol
- hook up column (Cobalt)—make sure there is no air in the system
- wash column with 5 column vol of ethanol, 5 col vol of water
- switch pumps to elution and wash buffers
- wash pumps

Adapted from Gelles Lab Protocol by Marc Ridilla and Linnea Lemma
Updated 2018-09-27 by Linnea Lemma

- wash column with 5 col vol elution buffer, 5 col vol lysis buffer. Do this at max speed possible within pressure limits.
- once there is a stable baseline in absorbance and conductivity – rezero

Cobalt Column—HiTrap TALON

- col vol is 5mL
- default flow rate is 5mL/min
- max flow rate is 20mL/min
- max pre-col pressure is 0.5MPa

Nickle Column—HisTrap HP

- col vol is 1mL
- default flow rate 1mL/min
- max flow rate 4.0mL/min
- max pre-column pressure 0.5MPa

Load the Protein

- remove the supernatant from centrifuge tubes
- while AKTA is in Manual load, use a syringe to load the clarified lysate into the injection loop (be careful to avoid bubbles!!!)
- AKTA setting for loading 5mL of protein
 - 1/5th of default flow rate
 - Injection valve—inject mode
 - 5mL timer
- after the protein is loaded: run 15 col vol--or until there is a stable absorbance--of lysis buffer through the system up to default flow rate
- put eppendorfs in fraction collector
- change to manual load mode to not include the loop
- then elute at 1mL/min for 5 col vol: COLLECT FRACTIONS

Determining where the protein is and how much you have

- Mini Bradford: Determine which wells are blue—those are the one with protein which we should run on a gel and then combine and keep.
- Run a gel (SDS page, 8 or 12% agarose)

Concentrate and Store the protein

- use filter concentrators and swinging bucket centrifuge (3000G)
- blank with the stuff that come through the filter (“filtrate”) and measure concentration using nano drop or other uv-vis spectrometer
- concentrate as desired
- store in 30% (final) sucrose or glycerol
- flash freeze in liquid nitrogen

Clean up Akta

- (1) 5x column volumes of elution buffer over column
- (2) Switch pumps to H2O and ethanol

Adapted from Gelles Lab Protocol by Marc Ridilla and Linnea Lemma
Updated 2018-09-27 by Linnea Lemma

- (3) Wash pumps
- (4) Change to inject mode to include loop
- (5) 5x col vol of H₂O
- (6) 5x col vol of ethanol