

Dogic Lab Protocol for Transforming and Expressing Proteins in *E. Coli*

Day 1: Transformation

- (1) Use Addgene (<http://www.addgene.org/>) to figure out the necessary antibiotic for the plasmid and how much to add
- (2) Make media for plates
 - 10g/L sodium chloride
 - 5g/L yeast extract
 - 10g/L tryptone
 - 15g/L agar for plates to make gel
 - each plate needs ~30mLProcedure:
 - Use a beaker and stir bar
 - add < (total volume) of water
 - weigh and add media components
 - add agar straight to bottle (with stir bar) for autoclaving
 - add well stirred solution
 - leave cap loose
 - add indicator tape
 - autoclave (liquid, 45min)
- (3) Add antibiotics
 - make stock antibiotic (at 1000x concentration) if needed
 - allow media to cool (~30min) to touch
 - don't want agar to harden
 - needs to be cool enough to not damage the antibiotic (50C)
 - add antibiotic at required concentration
 - pour into plates
 - allow agar to set (~1hr—not time sensitive)
- (4) Transform cells
 - get competent cells – put on ice (Rosetta 2 plysS)
 - get DNA – put on ice
 - in a sterile environment: add 1uL of DNA to each competent cell eppendorf (10uL)-- don't mix, pipette into the middle
 - leave on ice for ≥ 5 min
 - heat shock at 42.5C for 45sec
 - put on ice for “recovery” 2 min
 - in a sterile environment: add 300uL of antibiotic-free media
 - put in incubator (37C) for ≥ 1 hr
- (5) Spread cells on plates
 - label plates (date, DNA, cells used, Hi/Lo, initials)
 - in a sterile environment:
 - sterilize cell spreader using eth. and flame
 - add 50uL of starter onto Lo and the rest onto Hi for each DNA

allow to cool (antibiotics get damaged above ~50C)
spread cells well around plates
store upside down in incubator (37C) overnight

Day 2: Prepare Flasks, Begin starter growth

- (1) Check plates for growth
 - negative control has no growth
 - Hi has more colonies than Lo
 - seal one plate from each DNA (either Hi or Lo—want several colonies but dilute enough to pick a single one) and place in cold room to halt growth
- (2) Make media (2XYT)
 - 16g/L tryptone
 - 10g/L yeast
 - 5g/L NaCl
 - For BCCP plasmids, add 24mg/L (100uM) biotin
 - Procedure:
 - the largest beaker can hold 5L so concentrate as needed from there
 - use big beaker and large stir bar
 - make enough media for growth flasks and starter flasks
 - add foam stoppers, cover with foil & indicator tape
 - autoclave for 45min
- (3) Begin starter growth
 - allow media to cool
 - add antibiotics (100ug/mL Amp, 25ug/mL Chlor)
 - label flasks
 - scrape a single colony from each plate with pipette tip and place in each of the small flasks
 - place an unused pipette tip in negative control flask
 - put all flasks in incubator (37C, 200RPM)
 - make glycerol stock tubes (0.9mL 20% glycerol—sterilized by filter)

Day 3: Grow

- (1) Check that overnight cultures have growth (are cloudy)
- (2) Make glycerol stocks (0.9mL of 20% glycerol, 0.9mL of overnight—freeze at -80C for future use)
- (3) Inoculate broth with 1/200th volume of the overnight culture (2.5mL for 0.5L media)
- (4) Grow in air shaker (250rpm) at 37C and follow OD600
- (5) Induce expression at OD600=0.5-0.6.
 - remove flasks and place in cool water
 - add 1/1000th volume of 1M IPTG (final concentration is 1mM)
 - reset shaker thermostat to 18C
 - OPTIONAL: add 1/200th volume of 40mM Rifampicin 2hrs after induction (increases purity, decreases yield)
- (6) Shake overnight at 20C (~12-20hrs)

Day 4: Harvest

- (1) Pour cells into centrifuge bottles and spin at 5000G for 30min
- (2) Pour off most of the supernatant—save 50mL or whatever you want to resuspend in.
- (3) Resuspend pellet in 50mL of media and transfer to falcon tubes
- (4) Spin in swinging bucket centrifuge at 4C at maximum speed for 30 minutes.
- (5) Pour of supernatant. Freeze cells at -80C until ready to purify.
- (6) Rinse flasks with bleach. Rinse 3x volume with water. OPTIONAL: Autoclave dry for 60min.